



GEISTER®

The better way to operate™

**I N K A**™

**SURGICAL INSTRUMENTS**

Suppliers of Fine Quality Instruments

116 Cutler Road • Jandakot WA 6164 • AUSTRALIA  
customers@inkasurgical.com.au  
Phone: (08) 9417 9545 • Freecall: 1800 756 757 • Fax: (08) 9414 8906  
www.inkasurgical.com.au

GEISTER

Bayonet Micro Instruments  
for fine action in narrow and deep operative fields

## Rotatable shaft with flushing channel




41-4400 Micro Needle Holder, bayonet shaped, 26cm, rotatable, straight, jaw 0,5mm, Sapphire™ jaws, single action, WL 12cm

41-4402 Micro Needle Holder, bayonet shaped, 26cm, rotatable, straight, jaw 1,0mm, Sapphire™ jaw, single action, WL 12cm,




41-4401 Micro Needle Holder, bayonet shaped, 26cm, rotatable, curved, jaw 0,5mm, Sapphire™ jaw, single action, WL 12cm


41-4403 Micro Needle Holder, bayonet shaped, 26cm, rotatable, curved, jaw 1,0mm, Sapphire™ jaw, single action, WL 12cm




41-4410 Micro Cupped Forceps, bayonet shaped, 26cm, rotatable, straight, jaw Ø 2mm, double action, WL 12cm




41-4414 Micro Cupped Forceps, bayonet shaped, 26cm, rotatable, straight, jaw Ø 4mm, double action, WL 12cm




41-4411 Micro Cupped Forceps, bayonet shaped, 26cm, rotatable, 60° angled, jaw Ø 2mm, double action, WL 12cm,




41-4420 Micro Scissors, bayonet shaped, 26cm, rotatable, straight, sharp/sharp single action, WL 12cm,




41-4422 Micro Scissors, bayonet shaped, 26cm, rotatable, straight, blunt/blunt single action, WL 12cm,



41-4421 Micro Scissors, bayonet shaped, 26cm, rotatable, curved, sharp/sharp single action, WL 12cm,



41-4423 Micro Scissors, bayonet shaped, 26cm, rotatable, curved, blunt/blunt single action, WL 12cm,



41-4424 Micro Scissors, bayonet shaped, 26cm, rotatable, 45° angled, blunt/blunt single action, WL 12cm,



41-4426 Micro Scissors, bayonet shaped, 26cm, rotatable, 45° angled, sharp/sharp single action, WL 12cm,



41-4430 Micro Forceps, bayonet shaped, 26cm, rotatable, straight, jaw 0,5mm, single action, without ratchet, WL 12cm,



41-4431 Micro Forceps, bayonet shaped, 26cm, rotatable, curved, jaw 0,5mm, single action, without ratchet, WL 12cm,

Standard working length 12 cm, for 9 cm WL please add "S" behind article number. Additional working lengths on request.