

I N K A™

SURGICAL INSTRUMENTS

Suppliers of Fine Quality Instruments

116 Cutler Road • Jandakot WA 6164 • AUSTRALIA

customers@inkasurgical.com.au

Phone: (08) 9417 9545 • Freecall: 1800 756 757 • Fax: (08) 9414 8906
www.inkasurgical.com.au

GEISTER®

The better way to operate™



GEISTER

Basic Micro Set &
Nerve Repair Set



Basic Micro Set:

- 03-0055 Adventitia Scissors RH 15cm straight, short jaw
- 03-0056 Dissecting Scissors RH 15cm curved, short jaw
- 03-0111 Micro Needle Holder RH 15cm straight, without lock
- 03-0112 Micro Needle Holder RH 15cm curved, without lock
- 10-3226 Suture Tying Forceps 15cm straight
- 10-3227 Suture Tying Forceps 15cm curved
- 10-3215 Jeweller's Forceps 11cm straight (2 each recommended)
- 10-3200 Vessel Dilator Forceps 11cm 10° 0.2mm smooth jaws
- 75-0072 Clamp Applying Forceps Fig. 1-2 without catch
- 75-0325 Ven. Plain Single Clamps #4 1.5-3.5mm (pair)
- 75-0398 Ven. Plain Single Clamps #2 0.6-1.5mm (pair)
- 75-0400 Ven. Plain Single Clamps #3 1.0-2.25mm (pair)
- 77-0104 Micro Instrument Container 300x220x50mm
- 77-9017 Clamp Box 70x45mm

Our Basic Micro Set equips the surgeon with the basic instruments for most micro surgery procedures. The reduced Nerve Repair Set is suitable for Nerve and Tendon Repair.

Nerve Repair Set:

- 03-0060 Adventitia Scissors FH 15cm straight, long jaw
- 03-0112 Micro Needle Holder RH 15cm curved, without lock
- 10-3215 Jeweller's Forceps 11cm straight (2 each)
- 77-0101 Micro Instrument Container 220x165x50mm



© Geister 2010 · Printed in Germany · All rights reserved.



Geister Medizintechnik GmbH
Föhrenstrasse 2
D-78532 Tuttlingen/Germany
Phone +49 (0) 7461 – 966 240
Telefax +49 (0) 7461 – 966 2422
E-mail info@geister.com
Internet www.geister.com

Distributor